



ID-Centrifuge L · ID-Incubator L



Centrifuge and Incubator for ID-Cards and tubes





BIO-RAD ID-Centrifuge L 06-03-2014 10:28

00:00 0

BIO-RAD

ID-Centrifuge L











ID-Centrifuge L

Innovation

- Only one instrument for both ID-Cards and tubes
- Head for post transportation ID-Card centrifugation with adapted program (5 minutes at 85 g)
- Panel PC touch-screen
- Automatic head type detection and selection of the appropriate settings
- Multi-language software

Main Features

- Only one centrifuge for tubes and ID-Cards
- Automatic detection of the head type
- Automatic selection of the appropriate program
- Head for post transportation ID-Card centrifugation with adapted program (5 minutes at 85 g)
- 28 tube capacity
- Panel PC touch-screen with easy-to-use interface (icon based)
- 4 compatible centrifuge heads
- Multi-language software
- CE compliant according to IVD Directive 98/79/EC

Technical Data

Performance

RCF maximum	1,000 g (for tubes)
RCF	85 g (for ID-Cards)

Electrical data

Voltage	110-240 VAC
Frequency	50/60 Hz

Dimensions

Width	380 mm
Height	230 mm
Depth	500 mm
Weight	13 kg
(heads not included)	

Ordering Information

Catalog #	Description
-----------	-------------

ID-Centrifuge L

009201	ID-Centrifuge L , 1 unit
009410	Head 24 ID-Cards , 1 unit
009415	Head 28 tubes , 1 unit
009413	Head 12 tubes , 1 unit
009414	Head 2 racks , 1 unit



BIO-RAD

ID-Incubator L









ID-Incubator L

Innovation

- Only one instrument for both ID-Cards and tubes
- Incubation zones are independent and can be deactivated
- Panel PC touch-screen
- Multi-language software

Main Features

- Incubator for both tubes and ID-Cards
- Independent incubation areas (ID-Cards and/or tubes)
- Up to 36 ID-Card and 24 tube capacity
- Panel PC touch-screen with easy-to-use interface (icon based)
- Multi-language software
- Batch specific timers (3 for ID-Cards, 2 for tubes)
- Standard program: 15 minutes at 37°C
- Programmable incubation time
- CE compliant according to IVD Directive 98/79/EC

Ordering Information

Catalog # Description

ID-Incubator L

009203 **ID-Incubator L**, 1 unit

Technical Data

Electrical data

Voltage **110-240 VAC**
Frequency **50/60 Hz**

Dimensions

Width **380 mm**
Height **230 mm**
Depth **500 mm**
Weight **13 kg**
(heads not included)



Visit www.bio-rad.com/immunohematology for more information.



**Bio-Rad
Laboratories**

*Clinical
Diagnostics Group*

Website www.bio-rad.com/diagnostics **Australia** 61-2-9914-2800 **Austria** 43-1-877-8901 **Belgium** +32 (3)710-53-00 **Brazil** +55 (31)3689-6600
Canada 1-514-334-4372 **China** 86-21-61698500 **Czech Republic** 420-241-430-532 **Denmark** +45-4452-1000 **Finland** 358-9-804-22-00
France 33-1-47-95-60-00 **Germany** +49 (0)89-318-840 **Greece** 30-210-7774396 **Hong Kong** 852-2789-3300 **Hungary** +36-1-459-6100
India 1800-180-1224 **Israel** 972-3-9636050 **Italy** +39-02-216091 **Japan** 81-3-6361-7070 **Korea** 82-2-3473-4460 **Mexico** +52 (55)5488-7670
The Netherlands +31-318-540666 **New Zealand** 64-9-415-2280 **Norway** +47-23-38-41-30 **Poland** 48-22-3319999 **Portugal** 351-21-472-7700
Russia +7-495-721-1404 **Singapore** 65-6415-3170 **South Africa** 27-11-442-85-08 **Spain** 34-91-590-5200 **Sweden** 46-8-555-127-00
Switzerland +41 (0)26-674-55-05/06 **Taiwan** 886-2-2578-7189 **Thailand** 662-651-8311 **United Kingdom** +44 (0)20-8328-2000