




The power of syndromic testing

Fast and accurate answers
when it matters most

Sample to Insight

Syndromic solutions extend throughout your hospital system



Support rapid treatment decisions for patients

Timely and accurate results
+
Comprehensive panels with 21-24 targets
=
A shortened period of symptoms and reduced hospital stay



Empower confident diagnosis and therapy selection


Comprehensive results in about an hour = support rapid decisions on treatments

Ct values and amplification curves = potential insights into coinfection and disease severity



Improve antimicrobial stewardship at your hospital

Accurate, comprehensive results = target pathogen with initial antimicrobial therapy = reducing readmission and preventing hospital acquired infections



Provide more clinical insights with a streamlined workflow

Less than a minute hands-on time
+
Ct values and amplification curves
=
efficient delivery of valuable results, including insights on coinfection and sample quality



The next generation of syndromic insights

DiagCORE[®] is the next generation of syndromic insights you have been waiting for. Powered by proven QIAGEN sample and assay technologies, you can now confidently provide clinical insights to patients with the highest level of versatility.



Results in
about one
hour



Accurate and
comprehensive
results



Less than
1 minute
hands-on time



Ordering Information - Instrument

| Product | Contents | Product code |
|-----------------------------|---|--------------|
| DiagCORE Analyzer 1.0 | Instrument consists of 1 each of both the Operational Module and Analytical Module. | 100000 |
| DiagCORE Analytical Module | 1 each of module containing hardware and software for sample testing and analysis | 110001 |
| DiagCORE Operational Module | 1 each of module to enable interaction with the Analyzer | 110002 |

Ordering Information - Cartridge

| Product | Contents | Product code |
|------------------------------------|--|--------------|
| DiagCORE Respiratory Panel V2 | 6 individually packaged cartridges containing all reagents needed for sample preparation and multiplex RT-real time PCR plus internal control, including 6 transfer pipettes | 210002 |
| DiagCORE Gastrointestinal Panel V2 | 6 individually packaged cartridges containing all reagents needed for sample preparation and multiplex RT-real time PCR plus internal control, including 6 transfer pipettes | 210004 |

DiagCORE is intended for in vitro diagnostic use in Europe.

*QIAstat-Dx coming soon. Currently available as DiagCORE. Not available in countries.

Discover QIAstat-Dx – the next generation of syndromic insights

Visit QIAstat-Dx.com for more information

Trademarks: QIAGEN®, Sample to Insight®, QIAstat-Dx®, DiagCORE® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1116104 04/2019 PROM-13621-001 © QIAGEN 2019, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com

