



Looking for a cost effective solution to replace your old thermal cyclers?

## Learn more about the QIAamplifier® for end point cycling in human identity

Enjoy success at first try with high performance QIAGEN® Multiplex STR amplification and instrumentation


Discover Multiplex STR amplification with Quality Sensor for human identification and paternity testing



Fast cycling: 30 cycles in less than 45 min



Whisper quiet: low noise emission (45 dB max.)



Linear Gradient Function



Designed for robustness and reproducibility of results



# The QIAamplifier is a versatile system designed for performance

The QIAamplifier 96 is a high-performance thermal cycler with a 96-well block and a gradient function ideal for all your end-point PCR requirements.

- Easy-to-use: intuitive software & a 7-inch color touchscreen interface
- Fast ramping: heating and cooling rates of up to 4°C/sec
- Fast cycling: 30 cycles in less than 45 minutes
- High-performance smart lid: pressure control for highly reproducible results
- Whisper quiet: low noise emission (45 dB max.)
- Linear Gradient Tool: linear gradient from 20°C to 99°C (max. span 20°C and minimal increment of 0.1°C)

Contact your QIAGEN HID sales rep.



QIAGEN is committed to accelerate your research and provide your lab with high-quality PCR solutions meeting all your needs. Save time, money and peace of mind by choosing optimized AllTaq chemistries and a QIAamplifier cycler for fast, easy and powerful analysis of your PCR experiments.

Trademarks: QIAGEN®, Sample to Insight®, QIAamplifier® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

PROM-17193-001 1122845 10/2020 © 2020 QIAGEN, all rights reserved.

Ordering [www.qiagen.com/shop](http://www.qiagen.com/shop) | Technical Support [support.qiagen.com](mailto:support.qiagen.com) | Website [www.qiagen.com](http://www.qiagen.com)